REF G-24040, G-24041, G-24046, G-24047, G-24048, G-24049, G-24050, G-24051, G-24055, G-24056, G-24060, G-24061, G-24065, G-24070, G-24072, G-24310, G-24313, G-24315, G-24360, G-24362, G-24364, G-24366

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**ULTRASONIC TIPS
PRODUCT RELATED PROCESSING INSTRUCTIONS**

These reprocessing instructions for critical* medical devices (MD), are exclusively valid for the above listed items by GEUDER AG. All medical devices that are provided in a non-sterile way must be reprocessed before each application.

Fig. 1: example illustrations, ultrasonic tips

REF		to use with WRENCH		to use with GEUDER INFUSION SLEEVES														
		G-24110	G-24112	G-24073	G-24086	G-24101	G-24102	G-24103	G-24104	G-24106	G-24111	G-24311	G-24314	G-24316	G-24361	G-24363	G-24365	G-24367
G-24040		X	-	-	-	X	X	-	-	-	X	-	-	-	-	-	-	-
G-24041		X	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
G-24046		-	X	-	-	-	-	-	-	-	-	-	-	X	-	-	-	
G-24047		-	X	-	-	-	-	-	-	-	-	-	-	-	X	-	-	
G-24048		-	X	-	-	-	-	-	-	-	-	-	-	-	-	X	-	
G-24049		-	X	-	-	-	-	-	-	-	-	-	-	-	-	-	X	
G-24050		X	-	-	-	X	X	-	X	-	X	-	-	-	-	-	-	
G-24051		X	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
G-24055		X	-	-	-	X	X	-	X	-	X	-	-	-	-	-	-	
G-24056		X	-	-	-	X	X	-	X	-	X	-	-	-	-	-	-	
G-24060		X	-	-	-	X	X	-	X	-	X	-	-	-	-	-	-	

REF		to use with WRENCH		to use with GEUDER INFUSION SLEEVES														
		G-24110	G-24112	G-24073	G-24086	G-24101	G-24102	G-24103	G-24104	G-24106	G-24111	G-24311	G-24314	G-24316	G-24361	G-24363	G-24365	G-24367
G-24061		X	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
G-24065		X	-	-	-	-	-	X	-	X	-	-	-	-	-	-	-	-
G-24070		X	-	-	-	-	-	X	-	X	-	-	-	-	-	-	-	-
G-24072		-	X	X	X	-	-	-	-	-	-	-	-	-	-	-	-	-
G-24310		X	-	-	-	-	-	-	-	-	X	-	-	-	-	-	-	-
G-24313		X	-	-	-	-	-	-	-	-	-	X	-	-	-	-	-	-
G-24315		X	-	-	-	-	-	-	-	-	-	-	X	-	-	-	-	-
G-24360		-	X	-	-	-	-	-	-	-	-	-	-	X	-	-	-	-
G-24362		-	X	-	-	-	-	-	-	-	-	-	-	-	X	-	-	-
G-24364		-	X	-	-	-	-	-	-	-	-	-	-	-	-	X	-	-
G-24366		-	X	-	-	-	-	-	-	-	-	-	-	-	-	-	-	X

Tab. 1

ACCESSORIES

G-24110	WRENCH for reusable ultrasonic tips	
G-24112	WRENCH for CMP ultrasonic tip	
G-31021	CLEANING ADAPTER for ultrasonic tips	
ULTRASONIC HANDPIECES	example illustration	
GEUDER INFUSION SLEEVES	example illustrations	

Tab. 2, illustrations not drawn to scale

OVERVIEW

Cycle limit	Assembly/Disassembly	Accessories
100 x	YES	YES (see tab. 2 Accessories)
Pre-cleaning	Ultrasound	Automated processing
YES	YES	YES
Manual processing	Compressed air	Steam sterilization
NO	YES	YES

Tab. 3

AUXILIARY EQUIPMENT

First treatment	Soft tissue (e.g. compress), single-use syringe (5 times with 10 ml each)
Pre-cleaning	Soft plastic brush or steamer (steam jet, max. 2 bar), single use syringe (5 times with 10 ml each) or cleaning gun (water pressure gun, max. 2 bar) optionally: ultrasound bath 35-40 kHz
Drying	Medical compressed air (oil-free compressed air, max. 2 bar)
Checking pH value	pH test strips
Maintenance	Instrument oil apt for sterilization, lint-free wipes

Tab. 4

ASSEMBLY/DISASSEMBLY:

- Screw the ultrasonic tip onto the ultrasonic hand-piece and tighten it with the wrench for ultrasonic tips (Fig. 2).
- Slide the infusion sleeve over the ultrasonic tip and screw it onto the thread of the ultrasonic handpiece. Make sure not to damage the infusion sleeve (Fig. 3).
If you are using an infusion sleeve with a flexible core, make sure not to shift this core.
- For the correct positioning of the infusion sleeve refer to the product related processing instruction of the infusion sleeve (1102223).
- For the next steps please refer to the Instructions for Use for ultrasonic handpieces (1102217 IfU Ultrasonic Handpieces or 1102312 IfU Ultrasonic Handpiece ERGO)


Fig. 2


Fig. 3

FIRST TREATMENT AT PLACE OF USE

- Remove the infusion sleeve and the ultrasonic tip from the ultrasonic handpiece in reverse order as shown in fig. 2 to 3.
- The ultrasonic handpiece is reprocessed according to the respective instructions for use (1102217 or 1102312), the infusion sleeve according to its product related processing instruction (1102223).
- Remove macroscopically visible impurities (e.g. pharmaceuticals or corrosive solutions, etc.) during/ immediately after application with a compress soaked in sterile water.
- Rinse the ultrasonic tip with sterile water.
- Place dry ultrasonic tip into a suitable storage or transport system and close it.

PREPARATION BEFORE CLEANING

Pre-cleaning

- Rinse ultrasonic tip under running cold (< 40 °C) tap water.
- Remove visible, gross contamination using a soft plastic brush or a steamer.
- Attach the cleaning adapter as shown in fig. 4.
- Rinse ultrasonic tip with a cleaning gun with a suitable flushing adapter or a single-use syringe.
- Rinse ultrasonic tip against the flow direction.


Fig. 4

OR

Ultrasound cleaning as supporting pre-cleaning (optional)

Basically, we recommend pre-cleaning with ultrasound for all complex or very dirty medical devices (e.g. instruments with a lumen, slide switch, joint, etc.). Observe the instructions or exposure time and concentration of the producer of the chemical as well as the instructions of the manufacturer of the ultrasound bath.

- Introduce the ultrasonic tip and make sure that it is completely moistened:
 - Fill ultrasonic tip with cleaning agent.
- Ultrasound treatment at 35-40 kHz, wait until the exposure time is over.
- Rinse the ultrasonic tip afterwards with demineralised water for neutralisation.
- Check ultrasonic tip for unobstructed flow:
 - Rinse the ultrasonic tip with a cleaning gun, a steamer or a single-use syringe.
 - Rinse ultrasonic tip against the flow direction.

AUTOMATIC CLEANING AND DISINFECTION

As additional load only ophthalmic products are permitted.
Thermic disinfection A₀ value > 3000, 93 °C, 5 min; VARIO TD program.
The selection of suitable racks, charging trolleys for optimal rinsing and flushing of the medical devices should be done after consulting with the cleaning device manufacturer.

Process

- Connect the cleaning adapter to a suitable rinsing connector of the charging trolley.
 - Avoid rinsing shadows.
- After cleaning and disinfection the cleaning adapter must be removed for the following steps of reprocessing.

DRYING

Dry the ultrasonic tip with medical compressed air, if necessary.

Control of pH value

We recommend carrying out a batch-related routine control for alkaline residues.

- Apply the remaining humidity of the ultrasonic tip onto a pH test strip.
- Compare the pH test strip with the color scale:
 - pH must be neutral.

CONTROL AND MAINTENANCE

Visual control

- Check if the ultrasonic tip is clean:
 - In case of residues, the cleaning and disinfection process must be repeated.
- Check if the ultrasonic tip is damaged:
 - Remove corroded or damaged (bent) ultrasonic tips.
(Bent ultrasonic tips must not be straightened.)
- Check if the ultrasonic tip is complete.
- Check if the ultrasonic tip is dry.

Maintenance

- Care for all ultrasonic tips with a suitable instrument oil:
 - Let the ultrasonic tip cool down.
 - Apply the instrument oil on the thread
 - Remove excess instrument oil, if any, with a lint-free wipe.

PACKAGING

Use packaging materials and systems complying with the following standards:
DIN EN (ANSI AAMI) ISO 11607 und DIN EN 868-8.

STERILIZATION

Steam sterilization in fractionated prevacuum method:

SAL 10⁻⁶, at least 3 min at 134 °C, 3 bar.

STORAGE

Storage is at room temperature (max. 25 °C) in closed storage systems (dust- and light-protected also dry), protected from damage or mechanical influences, in a controlled environment.

- Damaged or open sterile barrier systems (SBS) are deemed non-sterile.
- The storage time depends on the type of packaging and must be validated accordingly.

TRANSPORT

Transport must not have a negative effect on the characteristics of the sterile goods.

Returns to Geuder AG: Please note that we only accept processed (cleaned, disinfected and sterilised) medical devices.

DISPOSAL

The laws and provisions of each country must be complied with. Medical devices suspected to contain CJD (vCJD): observe the recommendations of the Robert-Koch-Institute (www.rki.de).

National standards, regulations and/or restrictions must also be considered!

GENERAL INFORMATION FOR PROCESSING

Detailed information for processing are provided in the manufacturer information for the processing of medical devices from GEUDER AG.

The life cycle of reusable medical devices ends, either due to material changes, material wear, damages, as well as functional changes or loss of function of the medical product, or after reaching the maximum number of reprocessing cycles in case of a cycle limitation as stated in the product related reprocessing instructions for the product (see table overview).

DETAILS OF VALIDATION

Summary of processes validated by GEUDER AG:

The selection of the parameters refers to worst case conditions.

Manual pre-cleaning, cleaning and disinfection:

Pre-cleaning:

Using tap water (according TrinkwV [German Drinking Water Ordinance]), cleaning gun or syringe, instrument cleaning brush with soft nylon bristles (15 mm).

Cleaning:

- In the ultrasound bath, 0.5 %, 10 min/30-40 °C, 35 kHz, Bandelin RK 510 H; after ultrasound treatment the instruments are placed for 2 min in deionized water and rinsed and/or operated.

Process chemicals (manual):

- Cleaning: Dr. Weigert: neodisher MidiClean forte, deionized water.
- Disinfection: Bode Chemie: Korsorex plus, deionized water.

Automatic cleaning and disinfection:

Cleaning:

- Program: shortened VARIO TD program, 0.2 %, 55 °C/10 min, final rinsing 3 min/deionized water, Miele & Cie KG: Washer-Disinfector type G7836, charging trolley: E450/1 accessories for MIC instruments, sieve basket with lid.

Disinfection:

- Thermal high level disinfection A₀-value > 3000, 93 °C/5 min, Miele & Cie KG: Washer-Disinfector type G7836, charging trolley: E450/1 accessories for MIC instruments, sieve basket with lid.
- Thermal disinfection A₀-value = 600, 90 °C/1 min, Miele & Cie KG: Washer-Disinfector type G7836, charging trolley: E450/1 accessories for MIC instruments, sieve basket with lid.

Process chemicals (automatic):

- Cleaning: Dr. Weigert: neodisher MediClean forte.
- Final rinsing: deionized water.
- Disinfection: deionized water.

Steam sterilization:

- Program: half-cycle SAL10⁻⁶, fractionated pre-vacuum procedure 3 times, 132 °C/1.5 min, MMM GmbH: type Selectomat PL/669-2 CL, sterilization pouches, drying time 20 min, cooling time: outside of the autoclave 30 min at room temperature (max. 25 °C).
- The maximum sterilization temperature is 138 °C.

An independent, accredited test laboratory proved the basic suitability of the medical devices for an effective reprocessing.

The above mentioned instructions were validated by the manufacturer of medical devices to be suitable for the processing/reprocessing of a medical product to be reused. The processing/reprocessing person(s) have the responsibility, that the actually applied procedure with the utilized equipment, materials and staff in the processing facility achieves the required result. Therefore the verification and/or validation and a routine monitoring of the procedure are indispensable.

Actual version available at GEUDER AG.

Subject to changes.