


UNO COLORLINE TROCAR SYSTEM

FOR SUTURELESS VITRECTOMY



23G

25G



SINGLE-USE TROCAR SYSTEM TI ONE STEP 23G

Uno Colorline GREEN: 23 gauge / 0.6 mm
3 titanium ports for one-step technique, incl. infusion tubing
G-42311 1 set per box, sterile
G-46311 6 sets per box, sterile

SINGLE-USE TROCAR SYSTEM TI ONE STEP 25G

Uno Colorline BLUE: 25 gauge / 0.5 mm
3 titanium ports for one-step technique, incl. infusion tubing
G-42511 1 set per box, sterile
G-46511 6 sets per box, sterile

- Easy incision thanks to a completely new geometry of the blade
- Maximum tightness of the valves
- Silicone Tip-Instruments usable
- Special surface for a tighter fit in the sclera
- Trocar handle with marker function at 3 and 4 mm
- Highest possible mobility of the bulbus thanks to the "Flathead"-Design
- Minimized risk of port loss
- Self-sealing scleroderma

Geuder[®]
Precision made in Germany

GEUDER AG reserves the right to make changes to technical details in response to recent developments. GEUDER does not assume liability for the accuracy of each individual statement. Illustrations not drawn to scale.

GEUDER AG · Hertzstrasse 4 · 69126 Heidelberg · Germany
Phone: +49 6221 3066 · Fax: +49 6221 303122 · www.geuder.com